



Memorandum

Date: March 13, 2009

To: All PBME staff

From: Stan Grzybek on behalf of PBME Equipment Management Committee

Subject: Increased Scheduled Maintenance Interval for GE Model MAC 5000 and GE Model MAC 5500 EKG Machines

The EMC has decided to increase the Scheduled Maintenance interval for the device type Electrocardiograph Machine for all GE MAC 5000 and MAC 5500 models from 12 months to 24 months.

This change applies to BWH, DFCI, NSMC and MGH.

The rationale for the increased interval is attached.
Effective April 1, 2009 the next scheduled SM date will be updated to reflect the new interval.

Direct any issues to Stan Grzybek

Committee Approval dates

PBME EMC: 12/11/2008

BWH MEMC: Pending

MGH MEMC: January 12, 2009

This Memo and supporting documents can be found on

http://biomed.partners.org/main/EquipmentManagement/MAC5000_5500_IntervalIncrease.pdf

L:\DBEMAIN\Equipment\Equipment Management\PBME EMC Committee\Topics